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
Sample ID	Location	QC Status	QC Message	MTB Status	Drug Resistant Status
TB 724	20230325*F21 20	Pass		Positive	Rifampicin
TB 725	20230325*B20 20	Fail	Not all Control Ass	Weak Positive	Undetermined
TB 726	20230325*D20 20	Pass		Negative	Undetermined
TB 727	20230325*A23 20	Pass		Positive	Rifampicin/Rifabutin Isoniazid
TB 728	20230325*B19 20	Pass		Positive	Streptomycin

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- Iyer, A., Ndlovu, Z., Sharma, J., Mansoor, H., Bharati, M., Kolan, S., Morales, M., Das, M., Issakidis, P., Ferlazzo, G., Hirani, N., Joshi, A., Tipre, P., Sutar, N., & England, K. (2023). Operationalising targeted next-generation sequencing for routine diagnosis of drug-resistant TB. *Public Health Action*, 13(2), 43–49. <https://doi.org/10.5588/pha.22.0041>
- Liang, R., Li, J., Zhao, Y., Qi, H., Bao, S., Wang, F., Duan, H., & Huang, H. (2024). A comparative study of MassARRAY and GeneXpert assay in detecting rifampicin resistance in tuberculosis patients' clinical specimens. *Frontiers in Microbiology*, 15. <https://doi.org/10.3389/fmicb.2024.1287806>
- Low, J. L., Wu, M. L., Aziz, D. B., Laleu, B., & Dick, T. (2017). Screening of TB actives for activity against nontuberculous mycobacteria delivers high hit rates. *Frontiers in Microbiology*, 8(AUG). <https://doi.org/10.3389/fmicb.2017.01539>
- Mansoor, H., Hirani, N., Chavan, V., Das, M., Sharma, J., Bharati, M., Oswal, V., Iyer, A., Morales, M., Joshi, A., Ferlazzo, G., Isaakidis, P., Ndlovu, Z., & England, K. (2023). Clinical utility of target-based next-generation sequencing for drug-resistant TB. *International Journal of Tuberculosis and Lung Disease*, 27(1), 41–48. <https://doi.org/10.5588/ijtld.22.0138>
- Martinson, N. A., Hoffmann, C. J., & Chaisson, R. E. (2011). Epidemiology of tuberculosis and HIV: Recent advances in understanding and responses. *Proceedings of the American Thoracic Society*, 8(3), 288–293. <https://doi.org/10.1513/pats.201010-064WR>
- Pai, M., Nicol, M. P., & Boehme, C. C. (2016). *Tuberculosis Diagnostics: State of the Art and Future Directions*. <https://doi.org/10.1128/microbiolspec>
- S., S. (2012). Diagnostic Methods for Mycobacterium tuberculosis and Challenges in Its Detection in India. In *Understanding Tuberculosis - Global Experiences and Innovative Approaches to the Diagnosis*. InTech. <https://doi.org/10.5772/31714>
- Wu, X., Tan, G., Yang, J., Guo, Y., Huang, C., Sha, W., & Yu, F. (2022). Prediction of Mycobacterium tuberculosis drug resistance by nucleotide MALDI-TOF-MS. *International Journal of Infectious Diseases*, 121, 47–54. <https://doi.org/10.1016/j.ijid.2022.04.061>
- Yang, H., Li, A., Dang, L., Kang, T., Ren, F., Ma, J., Zhou, Y., Yang, Y., Lei, J., & Zhang, T. (2023a). A rapid, accurate, and low-cost method for detecting Mycobacterium tuberculosis and its drug-resistant genes in pulmonary tuberculosis: Applications of MassARRAY DNA mass spectrometry. *Frontiers in Microbiology*, 14. <https://doi.org/10.3389/fmicb.2023.1093745>
- Yang, H., Li, A., Dang, L., Kang, T., Ren, F., Ma, J., Zhou, Y., Yang, Y., Lei, J., & Zhang, T. (2023b). A rapid, accurate, and low-cost method for detecting Mycobacterium tuberculosis and its drug-resistant genes in pulmonary tuberculosis: Applications of MassARRAY DNA mass spectrometry. *Frontiers in Microbiology*, 14. <https://doi.org/10.3389/fmicb.2023.1093745>

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