


Lucigen[®]

Simplifying Genomics

epicentre[®]

MasterPure[™] Yeast RNA Purification Kit

Cat. No. MPY03100

1. Introduction

The MasterPure™ Yeast RNA Purification Kit provides all of the reagents needed to purify RNA from cell types including: *Candida*, *Saccharomyces*, *Schizosaccharomyces* and filamentous fungi. The kit utilizes a rapid desalting process¹ to remove contaminating macromolecules, avoiding toxic organic solvents, bead-beating, and spheroplasting.

2. Product Specifications

Storage: Store the Proteinase K and RNase-Free DNase I from the MasterPure Yeast RNA Purification Kit at –20°C in a freezer without a defrost cycle. The rest of the kit may be stored at room temperature for ease of use.

Storage Buffer: RNase-Free DNase I is supplied in a 50% glycerol solution containing 10 mM Tris-HCl (pH 7.5), 10 mM MgCl₂, and 10 mM CaCl₂. Proteinase K is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 0.1 mM EDTA, 10 mM CaCl₂, 0.1% Triton® X-100, and 1 mM dithiothreitol (DTT).

Quality Control: The MasterPure Yeast RNA Purification Kit is function-tested by purifying RNA from *Saccharomyces cerevisiae*. RNA quality and yield are assayed by agarose gel electrophoresis, spectrophotometry and use as a template for RT-PCR.

3. Kit Contents

| Desc. | Concentration | Quantity |
|--|---------------|----------|
| The MasterPure Yeast RNA Purification Kit is available in a 100-purification size: | | |
| Extraction Reagent for RNA | | 60 mL |
| MPC Protein Precipitation Reagent | | 50 mL |
| Proteinase K | @ 50 µg/µL | 200 µL |
| RNase-Free DNase I | @ 1 U/µL | 500 µL |
| 2X T & C Lysis Solution | | 20 mL |
| RiboGuard™ RNase Inhibitor | @ 40 U/µL | 100 µL |
| 10X DNase Buffer | 2 mL | |
| (330 mM Tris-HCl [pH 7.8], 660 mM potassium acetate, 100 mM magnesium acetate, and 5 mM DTT) | | |
| TE Buffer | | 7 mL |
| (10 mM Tris-HCl [pH 7.5], 1 mM EDTA) | | |

4. RNA Purification Protocol

A. RNA Purification

1. Dilute 1 µL of 50 µg/µL Proteinase K into 300 µL of Extraction Reagent for RNA for each sample. A premix may be prepared for multiple samples.
2. Pellet cells by centrifugation and discard the supernatant. The optimal number of cells varies with the species, but 1-1.5 mL of a mid-log culture gives good results for many types of yeast.
3. Vortex mix 10 seconds to disperse the cell pellet.

4. Add 300 μL of Extraction Reagent for RNA containing the Proteinase K and mix thoroughly by vortexing.
5. Incubate at 70°C for 10-15 minutes; vortex mix every 5 minutes. Shorter incubations of 5-10 minutes can yield nearly as much RNA.
6. Place the samples on ice for 3-5 minutes and add 175 μL of MPC Protein Precipitation Reagent to 300 μL of lysed sample (solution may become cloudy). Vortex vigorously for 10 seconds.
7. Pellet the debris by centrifugation for 10 minutes at 4°C at $\geq 10,000 \times g$ in a microcentrifuge.
8. Transfer the supernatant fluid to a clean microcentrifuge tube and discard the pellet.
9. Add 500 μL of isopropanol to the recovered supernatant fluid. Invert the tube 30-40 times.
10. Pellet the RNA by centrifugation at 4°C for 10 minutes at $\geq 10,000 \times g$ in a microcentrifuge.
11. Carefully pour off or aspirate the isopropanol without dislodging the RNA pellet. If removal of contaminating DNA is required, proceed with DNase I treatment in part B, otherwise, continue with this protocol. The MasterPure Yeast RNA Purification Kit extracts yeast DNA much less efficiently than RNA, and DNase I treatment may not be needed.
12. Rinse twice with 70% ethanol, being careful to not dislodge the pellet. Centrifuge briefly if the pellet is dislodged. Remove all of the residual ethanol with a pipet.
13. Resuspend the RNA in 35 μL of TE Buffer.
14. Quantitate RNA by its absorbance at 260 nm, and obtain an A_{260}/A_{280} ratio. Yields have been in the range of 25-50 μg of RNA per mL of $A_{600}=1.0$ cultures of *S. pombe*, and A_{260}/A_{280} ratios have been greater than 2.0. Alternatively, determine RNA yields for small samples by fluorimetry.

B. Removal of Contaminating DNA from RNA Preparations

1. Remove all of the residual isopropanol with a pipet.
2. Prepare 200 μL of DNase I solution for each sample. Add 20 μL of 10X DNase Buffer to 175 μL of deionized water, then add 5 μL of RNase-Free DNase I.
3. Completely resuspend the nucleic acid pellet in 200 μL of DNase I solution.
4. Incubate at 37°C for 10 minutes.
Note: Additional incubation (up to 30 min) may be necessary to remove all contaminating DNA.
5. Add 200 μL of 2X T and C Lysis Solution; vortex mix for 5 seconds.
6. Add 200 μL of MPC Protein Precipitation Reagent (solution may become cloudy). Vortex mix 10 seconds; place on ice 3-5 minutes.
7. Pellet the debris by centrifugation at 4°C for 10 minutes at $\geq 10,000 \times g$ in a microcentrifuge.
8. Transfer the supernatant containing the RNA into a clean microcentrifuge tube and discard the pellet.

9. Add 500 μ L of isopropanol to the supernatant. Invert the tube 30-40 times.
10. Pellet the purified RNA by centrifugation at 4°C for 10 minutes in a microcentrifuge at $\geq 10,000 \times g$.
11. Carefully pour off or aspirate the isopropanol without dislodging the RNA pellet.
12. Rinse twice with 70% ethanol, being careful to not dislodge the pellet. Centrifuge briefly if the pellet is dislodged. Remove all of the residual ethanol with a pipet.
13. Resuspend the RNA in 35 μ L of TE Buffer.
14. Add 1 μ L of RiboGuard™ RNase Inhibitor.

5. References

1. Miller, S.A. et al., (1988) *Nucl. Acids Res.* **16**, 1215.

Epicentre is a trademark of Illumina, Inc. and/or its affiliate(s) in the U.S. and other countries, and is used under license.

MasterPure and RiboGuard are trademarks of Lucigen.

Triton is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania. Vacutainer is a registered trademark of Becton Dickinson Corp., Rutherford, New Jersey.